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



Tel.: +49 (0) 22 24 / 818-140

Fax: +49 (0) 22 24 / 818-166

**Customer service /
Technical service
(Spare part orders for
components):**

JK-International GmbH,
Division JK-Global Service

Köhlershohner Straße
53578 Windhagen
GERMANY

Tel.: +49 (0) 22 24 / 818-863

Fax: +49 (0) 22 24 / 818-205

E-Mail: service@jk-globalservice.de

NOTE:

The illustrations and text apply to all units, even if only one is shown unless explicit differentiation is made between wellsystem MEDICAL / MEDICAL_PLUS and wellsystem RELAX / RELAX_PLUS.

NOTE:

The right to make technical modifications concerning the specifications and information in these installation instructions is reserved.

Reprinting and reproduction, in whole or in part, of this manual is only permitted with our previous written approval and with specification of the source.

Copyright

The copyright remains with JK-Holding GmbH.

The contents may not be partially or completely duplicated, distributed or utilised for competition purposes without authorisation or disclosed to other parties.

Table of contents

1	General.....	5
1.1	Symbol explanation.....	5
2	Safety information for dismantling and assembly	6
3	Information for service staff.....	8
3.1	Transport and transport damage.....	8
3.2	Installation location.....	9
3.3	Electrical connections	11
3.4	Water supply	12
3.5	Test run/Commissioning	12
3.6	Maintenance	12
3.7	Disconnect the device from the power supply and ensure that it cannot be switched on again.....	13
3.8	Disposal	14
3.8.1	Environmental regulations – disposal of batteries.....	14
3.8.2	Operating water disposal	14
3.8.3	Packaging	14
3.8.4	Disposal of old devices	14
3.8.5	Registration numbers	14
4	Assembly and disassembly instructions.....	15
4.1	Tools and equipment.....	15
4.2	Scope of delivery.....	15
4.2.1	Supplied accessories	15
4.2.2	Optional accessories.....	16
4.3	Packaging units and module allocation	16
5	Disassembly	17
5.1	Removing the panels and the frame	17
5.2	Removing the control unit	21
5.3	Removing the drive motors	24
5.4	Removing the pump	25
5.5	Removing the well.....	25
5.6	Disassembling the frame.....	28
6	Assembly	29
6.1	Assembling the frame	30
6.2	Fitting the well	31
6.3	Installing and connecting up the pump.....	33
6.4	Fitting the control box.....	34
6.5	Aligning the massage unit (empty well).....	36

6.6	Check whether the nozzle carriages move smoothly	37
6.7	Filling the well.....	43
6.8	Fitting the drive motors.....	44
6.9	Connecting up the cooling system.....	45
6.10	Assembling the frame.....	46
6.11	Fitting the handle.....	47
6.12	Fitting the control panel (wellsystem RELAX / RELAX_PLUS).....	48
6.13	Fitting the control panel (wellsystem MEDICAL / MEDICAL_PLUS)	49
6.14	Fitting the panels.....	51
7	Technical data	54
7.1	Dimensions wellsystem RELAX	54
7.2	Dimensions wellsystem RELAX_PLUS	55
7.3	Dimensions wellsystem MEDICAL	56
7.4	Dimensions wellsystem MEDICAL_PLUS.....	57
7.5	Power & connection ratings.....	58
7.6	Weights	58
7.7	Water/cooling conditions	59
7.8	Coin devices.....	59
8	Appendix.....	60
8.1	Warning sticker.....	60
8.2	Controller overview.....	61
8.3	Wellsystem MEDICAL / MEDICAL_PLUS: EMC information	63
8.3.1	Table 201: Guidance and manufacturer's declaration – electromagnetic emissions	63
8.3.2	Table 202: Guidance and manufacturer's declaration – electromagnetic immunity.....	64
8.3.3	Table 204: Guidance and manufacturer's declaration – electromagnetic immunity.....	65
8.3.4	Table 206: Recommended separation distances between portable and mobile RF communications equipment and the massage unit	66
9	Index.....	67

1 General

1.1 Symbol explanation

The following types of safety notices are employed in these installation instructions:

DANGER!

Type and source of hazard

This safety notice indicates the existence of a direct danger to body and life.

DANGER!

Type and source of hazard

This safety notice warns of dangers to body and life which are caused by electricity.

WARNING!

Type and source of hazard

This safety notice warns of machine, material or environmental damage.

NOTE:

This symbol does not identify any safety notices, but provides additional information to better understand the processes.

2 Safety information for dismantling and assembly

DANGER!

Electricity throughout entire device!

Assembly and connection to the mains electrical supply must comply with the national regulations.

wellssystem MEDICAL / MEDICAL_PLUS is subject to special preventive measures with regard to electromagnetic compatibility and must be installed and commissioned in compliance with the EMC instructions in the Appendix; see page 63.

Assembly, installation, upgrading and servicing the device is only to be carried out by personnel who have been specially trained by JK-Global Service.

- The device must only be set up and used in closed rooms.
- Safety equipment (e.g. switches) and safety instructions must not be removed or deactivated as this might impair a safe operation of the device!
- You must always abide by the danger and safety instructions warnings attached to the device, even during assembly and disassembly (see page 60).
- The unit is only to be handed over to the operator in a flawless condition!
- The massage unit must not be used if it has not been filled!
- A water cooling system must be connected to the device if it is to be used for more than 1 hour per day at maximum pressure. The pressure hoses must not be removed.
- The device may only be operated with the pressure hoses included in the scope of delivery.
- A water shut-off valve must be used.
- The assembly is only to be carried out in the sequence stipulated in these assembly instructions!
- wellssystem MEDICAL / MEDICAL_PLUS:
The unit must be disconnected from the power supply before any assembly and disassembly work is carried out on the electrical units (see page 13). Exceptions are permissible only for function tests!
- wellssystem RELAX / RELAX_PLUS:
Before beginning work, pull out the mains plug or disconnect all live lines.

DANGER!

Electricity throughout entire device!

- High frequency mobile communication systems can impair the functioning of the device.
- Heat therapy equipment (microwave therapy units, etc.) installed in the surrounding areas might cause electromagnetic interference. Please contact our customer service department if this happens (see page 2).


WARNING!

Operating faults!

The device might be damaged.

- The device should not be installed and operated on the transport pallet. If installed on the transport pallet, optimal alignment and adjustment is not possible. Breakdown and malfunctions are unavoidable!

14250 / 0

3 Information for service staff


3.1 Transport and transport damage

WARNING!

The equipment might be damaged if the transport safety devices are missing!

The transport safety device (1) must be fitted to the pump if the massage unit has to be transported with the pump installed.

12990 / 0

Please inform the owner if you discover any transport damage when unpacking and installing the unit.

Please take note of the following information, so that any damage can be taken care of quickly and without any complications arising:

- Any recognisable damage must be reported in writing within a period of 4 working days – this includes the delivery date and being received by the transport company – to the transport company that delivered the shipment.
- The resulting damage must be made good by the transport company that delivered the shipment as, according to the general terms and conditions of the transport company, the transport risk is assumed by the person who ordered the equipment. Naturally, we will fully support you in asserting your justifiable claim.

3.2 Installation location

WARNING!

Salty air! Contaminated air!

Device damage due to corrosion of the housing and electrical components.

- Do not install the device at swimming pools.
- Do not operate the device in locations with aggressive environmental conditions (e.g. air containing chlorine)

Damage to device possible!

Where there are great differences in temperature between the transport route and the installation site the device may not be put into operation immediately after being set up.

- Wait at least 2 hours before connecting the device to the power supply.

Risk of overheating due to air flow problems!


Operating faults possible.

- Do not change, obstruct or block the air inflow and outflow to the device or make any unauthorized changes to the device. The manufacturer shall bear no liability for damage or injuries resulting from such actions.
- Observe the ventilation system data.
- Do not set up and operate the device on the transport pallet!

Permissible Floor / Ceiling Loads

The massage unit has been designed to be installed on a flat and level surface.

When installing the unit, you should always make sure that the floors and ceilings in commercially used rooms are designed to support a load of 5000 N/m².


NOTE:

The load-bearing capacity of wooden beam ceilings must be proven in individual cases.

If the actual maximum load exceeds this value, the operator must provide separate documentary evidence in conformity with DIN EN 1991-1-1:2010-12 for the use of these rooms.

Example for ceiling load (based on normal room or cabin dimensions, 1 device including 1 user, 1 member of staff and small furnishings):

1 device up to 800 kg: Ceiling load approx. 2150 N/m².

The massage unit can also be installed on a plywood plate in order to better distribute the adjustable feet pressure over the ground.

L = 2150 mm

W = 1000 mm

Please contact the customer service department without fail before putting the device into operation - see page 2.

Environmental conditions


- Optimum ambient temperature: Between 15 °C and 32 °C
- Storage temperature: Between -15 °C and +60 °C (the well, pump and the hoses must not contain any water)
- Relative humidity: Between 30 % and 75 %
- Air pressure: Between 700 hPa and 1060 hPa

The device is not to be used outdoors.

Adaptation of the unit is required for operating altitudes greater than 2000 m above sea level in order to guarantee trouble-free operation.

The device is not suitable for operation in mobile facilities (ships, buses, trains). Adaptation of the unit is required in order to guarantee trouble-free operation when used in mobile facilities.

3.3 Electrical connections

- Only an authorised specialist company is permitted to make the on-site electrical connections (equipment connection socket)!
- You must abide by the regulations.
- We recommend installing a selective RCD - circuit breaker (Residual Current Device) 30 mA, type B, class C. The RCD - circuit breakers must be selective or short time delay.
- Only one massage unit should be connected to each RCD.
- wellsystem MEDICAL / MEDICAL_PLUS: A maximum mains impedance of 0.14 ohms must be guaranteed for pulse massage.
- The electrical installation must have a freely accessible, all-pole disconnecting device (Master switch with clearly labelled On / Off settings, e.g. '0' and 'I' or additional indicator light) fitted externally on the site side.
- The unit assignment can be seen from the inscription.
- wellsystem RELAX / RELAX_PLUS:
The socket which is required for the electrical connection must be freely accessible or fitted with an easily accessible all-pole isolating device (master switch), complying with overvoltage category III. This means that each pole must have a contact opening width complying with the conditions of overvoltage category III for full isolation.
- wellsystem MEDICAL / MEDICAL_PLUS: The equipment must be permanently connected to the public power supply via an all-pole main switch. One power supply circuit (power supply lead) may only power one device.
- The power cable cross-section must be chosen in accordance with the national regulations (VDE 0100 in Germany) e.g. power cable lengths and the connected rating used must be taken into consideration.
- A supply line is required for each unit.
- Data and/or control lines are to be laid with a minimum distance of at least 10 cm from the mains cables.
- The mains connection cable must be laid before the unit is installed.

3.4 Water supply

The massage unit must be filled with normal mains water without additives. Always comply with local regulations when connecting to the public water supply mains.

Water supply and cooling conditions

- Permissible water pressure (water inlet): 0.8 MPa (8 bar)
- Maximum cooling water temperature: 16 °C

The device may only be operated with the pressure hose included in the scope of delivery.

- Pressure hose 90906-.. (10 m, 3/8");
Inlet: 1.5 MPa (15 bar) operating pressure
- Hose 50909-.. (10 m, 1/2"); outlet

3.5 Test run/Commissioning

- Has it been filled with sufficient water?
- If daily operation will be more than 1 hour at maximum pressure: has a cooling system been attached?
- Is the heating connected up?
- Are you sure that there are no leaks coming from the tank or pump?

If all questions were answered with 'Yes', then you can start the test run or operate the unit.

- Switching the unit on – see operating instructions.

3.6 Maintenance


Remount any protective devices (e.g. safety switch) after completing work.

WARNING!

Use only genuine spare parts of the same type! The use of other parts voids the CE conformity! Any liability is excluded in the case of damage or injury proven to be the result of the use of non-genuine spare parts.

3.7 Disconnect the device from the power supply and ensure that it cannot be switched on again

DANGER!

Electricity throughout entire device!

If work is to be performed on the devices, disconnect them from the power supply. This means that all live cabling must be switched off.

wellsystem MEDICAL / MEDICAL_PLUS:

Disconnect the sunbed from the power supply when performing maintenance work requiring opening of the device.

wellsystem RELAX / RELAX_PLUS:

Before beginning work, pull out the mains plug or disconnect all live lines.


Switching off the device alone is insufficient, as certain points may still be live. Therefore, switch off all fuses before performing any work and – if possible – remove them.

Personal danger due to electric shock!

Severe accidents can occur if unintentionally switched on again.

Immediately after disconnecting from the power supply, ensure all switches or fuses used to disconnect the power supply can not be switched on again.

- Lock the fuse box with a padlock.

In case of circuit breakers which cannot be unscrewed, a strip of adhesive tape may be applied over the actuation lever with the wording 'Do not switch, danger!' (1 + 2).

01511 / 1

A prohibition sign with the wording: 'Work in progress!'

- 'Place:
- 'Sign may only be removed by:

must always be securely displayed immediately.

DANGER!

Personal danger due to electric shock!

Prohibition signs must not be fixed to or come into contact with live parts.

3.8 Disposal

3.8.1 Environmental regulations – disposal of batteries

Batteries contain heavy metal compounds.


Within the European Union the national implementation of the Waste Framework Directive 2008/98/EC applies. According to the national waste disposal laws and in accordance with the community waste regulations, proof must be provided of the proper disposal of batteries.

3.8.2 Operating water disposal

The massage unit must be filled with normal mains water without additives. The water can be drained off into the public drainage system after use.

3.8.3 Packaging

All packaging consists of 100 % recyclable materials. Packaging brought into circulation by the JK Corporate Group that is no longer required can be returned to the JK Corporate Group. Your partner agency or dealer will be happy to advise.


3.8.4 Disposal of old devices

The device has been manufactured using recyclable materials. The JK Corporate Group will provide you with information on the content or potential hazards of the materials used.

In accordance with the 2012/19/EC directive, the manufacturer is obliged to take back and dispose of certain electrical and electronic components.

The components and devices are labelled with the following symbol.

The device will, if requested, be properly disposed of by the JK Corporate Group. This service is available free of charge. Your partner agency or dealer will be happy to advise.

3.8.5 Registration numbers

Companies within the JK Corporate Group are registered as manufacturers in Germany and accept all obligations in respect of the German Electrical and Electronic Equipment Act.

Registry number JK-Products GmbH (Devices): WEEE-DE 62655951

4 Assembly and disassembly instructions

At least two fitters are required for the disassembly and assembly.

4.1 Tools and equipment

DANGER!

Infections can be spread by skin contact!

- Always wear safety gloves when working on a massage unit that has already been commissioned, as it is always possible that you will come into contact with the water that is being used (e.g. when changing the drive belt or the rubber mat).
 - You must use a disinfectant afterwards to clean your hands, e.g. Skinman Soft SMS 50.
-

The following tools are required for disassembly and assembly:

- No. 2 Philips screwdriver
- No. 3 Philips screwdriver
- 8 to 19 mm double open-ended wrenches
- 4, 5 and 6 mm Allen keys
- Water pump pliers or 46 and 55 mm open-ended wrenches
- Torque Wrench (15720612) 8 mm (nut)
- Ratchet wrench with PHZ-bit 2 and size 7, 8, 10 and 13 mm box spanners
- Battery-powered screwdriver
- Size 4 (5/32") flat-bladed screwdriver
- Spirit level (at least 800 mm long)
- 1/2" water hose with Gardena hose-piece
- Industrial safety gloves

4.2 Scope of delivery


4.2.1 Supplied accessories

- Technical documentation (folder with instructions, brochures for error codes and additional documentation)
- Service card 12668-..
- Headrest 801091-..
- Pressure hose 90906-.. (10 m, 3/8");
Inlet: 1.5 MPa (15 bar) operating pressure
- Hose 50909-.. (10 m, 1/2"); outlet

4.2.2 Optional accessories

- Chip card set (default massage programs)
- Chip card set (blank chip cards = PROGRAMMABLE)
- Extension cushion 34528100

4.3 Packaging units and module allocation

12991 / 3

5 Disassembly

This section describes the steps for disassembly in detail.

DANGER!

Electricity throughout the entire device!

Assembly and connection to the mains electrical supply must comply with the national regulations.

wellssystem MEDICAL / MEDICAL_PLUS:


Disconnect the sunbed from the power supply when performing maintenance work requiring opening of the device.

wellssystem RELAX / RELAX_PLUS:

Before beginning work, pull out the mains plug or disconnect all live lines.

Assembly, installation, upgrading or servicing the equipment is only to be carried out by personnel specially trained by JK-Global Service.

5.1 Removing the panels and the frame


12992 / 1

- Remove the packing.
- Release the front panel from the holder.

Tool: No. 2 Philips screwdriver

- Remove the control panel, handle, headrest pad and hoses* from the unit and place to one side.

* 2 pressure hoses (90906-.. and 50909-..)


13374 / 0

- Undo the hose clip and pull the hose out.
- Remove the panel on the rear from the holder and place to one side with the hose.

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

Removing the frame panels


12993 / 0

- Unscrew the screws.

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

- Remove the panels.


Removing the panels

12994 / 0

- Unscrew the 6 x screws at the head end and remove the panel.
- Unscrew the 4 x screws at the feet end. Now slide the panel out to the side.


Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

12995 / 0

- Undo the 2 x screw sets.


Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

13061 / 0

- Slide the latched-in panels out along the rails.

Removing the transport safety device


12996 / 0

- Release the transport safety device. Unscrew the screws to do so.
- Unscrew the screws at the head and feet ends afterwards and then remove the plate.

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2 and 4 mm Allen key


- Retain the transport safety device.

5.2 Removing the control unit


12997 / 0

- Unplug the ST17 and ST6 connectors.

12998 / 0

- Unplug the ST18 and ST20 connectors.


12999 / 0

- Unscrew the screws.


Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

- Remove the cover from the control box.


13000 / 0

- Unplug the X018-1 and X019-1 connectors.
- Push the cable through the opening in the control box.

13001 / 0

- Unplug the X01-1, X02-1 and X03-1 connectors.
- Push the cables through the openings in the control box.


13002 / 0

- Unscrew the screw [1] in the centre first.
- Undo the upper screws [2] afterwards.

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2


- Unhook the control box.

5.3 Removing the drive motors

13003 / 0


- Unplug the connectors.
- Undo the cable ties.

13004 / 0

- Unscrew the nuts.
- Tool: Ratchet-wrench with 10 mm hexagonal socket wrench or size 10 mm open-ended wrench
- Remove the drive motors.


5.4 Removing the pump

- Undo the hose clips at the pump and then pull the hoses out.
- Unscrew the screws.


Tool: Ratchet wrench with 8 mm hexagonal socket wrench and 4 mm Allen key

5.5 Removing the well


- Undo the hose clip and then pull the hose out.

Tool: Ratchet wrench with 8 mm hexagonal socket wrench or 8 mm open-ended wrench

13053 / 0


- Unplug the connectors.

13054 / 0


- Unscrew the screws.
Tool: 4 mm Allen key
- Remove the plate holding the drive motor.

Undo the heater plug connections

13055 / 0


- Unplug the connectors.

13006 / 0

- Unscrew the screws.
- Tool: Ratchet-wrench with 4 mm Allen key
- Remove the well.
 - Lift the frame off the transport pallet.
 - Dispose of the transport pallet correctly, see Page 14.

5.6 Disassembling the frame


- Unscrew the screws.

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

- Disassemble the frame.

6 Assembly

Assembly module overview


13064 / 3

6.1	Assembling the frame.....	30
6.2	Fitting the well	31
6.3	Installing and connecting up the pump	33
6.4	Fitting the control box	34
6.5	Aligning the massage unit (empty well).....	36
6.6	Check whether the nozzle carriages move smoothly	37
6.7	Filling the well	43
6.8	Fitting the drive motors	44
6.9	Connecting up the cooling system.....	45
6.10	Assembling the frame.....	46
6.11	Fitting the handle	47
6.12	Fitting the control panel (wellsystem RELAX / RELAX_PLUS)	48
6.13	Fitting the control panel (wellsystem MEDICAL / MEDICAL_PLUS).....	49
6.14	Fitting the panels	51

This section describes the steps for assembly in detail.


6.1 Assembling the frame

- Optional: Place a plywood plate underneath the unit, see Page 9.
- Assembling the frame.
- Fit and tighten up the 16 x screws (4.2 x 9.5).

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

6.2 Fitting the well


13009 / 0

- Place the well on top.
- Fit the 24 x screws (M6 x 20) and tighten up.

Tool: Ratchet-wrench with 4 mm Allen key


Plug the heater plug connections in

13056 / 0

- Insert the plug.


Fitting drive motor 1 (well)

13057 / 0

- Fit the plate holding the drive motor in the unit.
- Fit the 4 x screws (M6 x 16) and tighten up.


Tool: 4 mm Allen key

13058 / 0

Insert the plugs.

- 2-pin BU13
- 4-pin BU23


13059 / 0

- Fit the hose and tighten the hose clip.

Tool: Ratchet-wrench with 8 mm hexagonal socket wrench or 8 mm open-ended wrench

6.3 Installing and connecting up the pump


17223 / 0

- Fit the pump.
- Fit the 2 x screws (M6 x 20), spring rings and washers and then tighten up.
- Fit the hoses and tighten the hose clips with 10 Nm.

Tool: Ratchet wrench with 8 mm hexagonal socket wrench, 4 mm Allen key and torque wrench (15720612)


6.4 Fitting the control box

13013 / 0


- Hook the control box on the upper screws.
- Tighten the screws.
- Fit the centre screw (4.2 x 9.5) and tighten up.

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2


13014 / 0

- Push the cables through the openings in the control box.
- Insert plugs X01-1, X02-1 and X03-1.


13015 / 0

- Push the cable through the opening in the control box.
- Insert plugs X018-1 and X019-1.

13017 / 0


- Plug in the ST18 and ST20 connectors.

13018 / 0

- Plug in the ST17 and ST6 connectors.

6.5 Aligning the massage unit (empty well)


13020 / 0

- Using a spirit level, check if the device is horizontal.
- In event the device is not standing horizontally, turn the screws on the feet until the device has been aligned correctly.


Tool: 17 mm open-ended wrench (locking nut)

6.6 Check whether the nozzle carriages move smoothly


- Loosen the screws and remove the supports.
- Remove the aluminium strips and fold the rubber blanket backwards.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench or 10 mm open-ended wrench

- Unscrew the screws.


Tool: Ratchet-wrench with 10 mm hexagonal socket wrench

13243 / 0


- Check whether the nozzle carriage moves smoothly. Move the nozzle carriage back and forth.

If the nozzle carriage moves back and forth without any problems, proceed with assembly, see Assembly on page 41. If the nozzle carriage moves with difficulty then perform the following steps.


13244 / 0

- Push the nozzle carriage to the foot end.


- Loosen the screws 2-3 mm.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench


- Retighten the screws.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench

13247 / 0


- Push the nozzle carriage to the head end.

13248 / 0

- Loosen the screws 2-3 mm.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench


13249 / 1

- Retighten the screws.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench


Assembly

13250 / 0

- Position the nozzle carriage (1) above the lower support.
- Insert and tighten the screws.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench


- Fold the rubber blanket forwards.
- Insert the aluminium strips.
- Mount 8 supports and attach with screws.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench or 10 mm open-ended wrench

6.7 Filling the well

The well must be vented approximately 3 days after filling (see the operating instructions for this device).

- Connect up the water hose and open the water shut-off valve.

WARNING!


Equipment damage can be caused through using the wrong filling quantity!

- Take note of the maximum filling quantity.

- Close the water shut-off valve and disconnect the water hose.
- Attach the remaining 4 supports.
- Retighten the screws.

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench or 10 mm open-ended wrench


6.8 Fitting the drive motors

13011 / 0

- Fit the drive motors.
- Tighten up the 8 x nuts (M6).

Tool: Ratchet-wrench with 10 mm hexagonal socket wrench or 10 mm open-ended wrench

13012 / 0

- Use the cable ties to secure the cable in place.

Insert the plugs:


Drive motor 2

- 2-pin BU14
- 4-pin BU24

Drive motor 3

- 2-pin BU15
- 4-pin BU25

6.9 Connecting up the cooling system


WARNING!


- The device may only be operated with the pressure hoses included in the scope of delivery.

- Attach the hoses.
- Set the post-cooling time (see operating instructions of the device).

6.10 Assembling the frame


- Slide the panels in from above onto the holders and align them to the screw holes.
- Push the panels down until you hear them latch into the groove.

- Fit the 4 screws (4.2 x 13) and tighten them up.
- Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

6.11 Fitting the handle


13024 / 1

- Position the handle and use the 4 x screws (M6 x 20) to secure it in place.

Tool: 4 mm Allen key

6.12 Fitting the control panel (wellsystem RELAX / RELAX_PLUS)


13050 / 0

- Position the control panel and use the 2 x screws (M6 x 20) to secure it in place.

Tool: 4 mm Allen key

- Push the cable through the opening in the control box.
- Insert the X016-1 plug.
- Now run a test, see Page 12.


6.13 Fitting the control panel (wellsystem MEDICAL / MEDICAL_PLUS)

- Position the control panel and use the 2 x screws (M6 x 20) to secure it in place.

Tool: 4 mm Allen key

- Push the cable through the opening in the control box.
- Insert the X016-1 plug.
- Now run a test, see Page 12.


13016 / 0

- Use the straps to hold the cover plate in the guides of the control box.
- Hold the top of the cover plate and tighten up the 2 x screws (M6 x 20).

Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2


6.14 Fitting the panels

13025 / 1

- Fit the panel at the head end and insert and tighten up the 6 x screws (M6 x 16).
- Slide the panel in sideways in the guides at the feet end. Insert and tighten up the 4 x screws (M6 x 16).


Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

13026 / 1

- Fit the panels.
- Fit the 2 x screws (M6 x 16) plus locking washers and tighten up.


Tool: No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

13375 / 0

- Attach the panel to the holder on the rear side.
- Attach the hose to the pump and tighten the hose clip.


Tool: Size 4 (5/32") flat-bladed screwdriver, No. 2 Philips screwdriver or battery-powered screwdriver with PHZ-bit 2

13027 / 2

- Attach the front panel to the holder.

Tool: Size 4 (5/32") flat-bladed screwdriver



13028 / 1

- Place the headrest pad in position.
- Test run/Commissioning, see page 12.

7 Technical data

7.1 Dimensions wellsystem RELAX

14351 / 1

A1 = 150 mm

A2 = 100 mm

BK = 2450 mm



TK = 1900 mm

L = 2150 mm

B = 1060 mm

H = 570 mm



7.2 Dimensions wellsystem RELAX_PLUS

13030 / 1

A1 =	70 mm
A2 =	70 mm
BK =	2450 mm
TK =	1900 mm
L =	2150 mm
B =	1060 mm
H =	570 mm



7.3 Dimensions wellsystem MEDICAL

14140 / 1

A1 =	150 mm
A2 =	100 mm
BK =	2450 mm
TK =	1900 mm
L =	2150 mm
B =	1100 mm
H =	570 mm

7.4 Dimensions wellsystem MEDICAL_PLUS

13029 / 1

A1 =	70 mm
A2 =	70 mm
BK =	2450 mm
TK =	1900 mm
L =	2150 mm
B =	1100 mm
H =	570 mm

7.5 Power & connection ratings

Device type:	wellssystem RELAX
Rated power consumption:	2600 W
Rated voltage:	230 - 240V AC
Type of connection:	Shock-proof plug
Rated frequency:	50 Hz
Fuses rating:	External: 16A

Device type:	wellssystem RELAX_PLUS
Rated power consumption:	2600 W
Rated voltage:	230 - 240V AC
Type of connection:	Shock-proof plug
Rated frequency:	50 Hz
Fuses rating:	External: 16A

Device type:	wellssystem MEDICAL
Rated power consumption:	2800 W
Rated voltage:	230 - 240V AC
Type of connection:	Fixed connection (3.00 m)
Rated frequency:	50 Hz
Fuses rating:	External: 16A

Device type:	wellssystem MEDICAL_PLUS
Rated power consumption:	2800 W
Rated voltage:	230 - 240V AC
Type of connection:	Fixed connection (3.00 m)
Rated frequency:	50 Hz
Fuses rating:	External: 16A

7.6 Weights

Empty weight (unfilled):	approx. 168 kg
Total weight (filled:)	approx. 472 kg

7.7 Water/cooling conditions

The massage unit must be filled with normal mains water without additives.
The water can be drained off into the public drainage system after use.

Permissible water pressure (water inlet):	0.8 MPa (8 bar)
Maximum coolant temperature:	16 °C > setpoint temperature and set time.


The device may only be operated with the high-pressure hoses included in the scope of delivery.

7.8 Coin devices

You can obtain information regarding the coin devices from the customer service department – see page 2.


8 Appendix

8.1 Warning sticker

13065 / 2

8.2 Controller overview


13060 / 0

Description	Information
BU21	JK coin device / external coin device connection
SD card	SD card insertion slot (software update)
ST1	LAN connection (for service personnel only)
ST6	Frequency converter fan
ST17	Frequency converter on the water pump
ST18	Heater
ST20	Limit switch, solenoid
X01-1	Drive motor (1) beneath the well (nozzle trolley)
X02-1	Drive motor (2) on feet side (front nozzle)
X03-1	Drive motor (3) on feet side (rear nozzle)
X04-1	Circuit board and drive motor power supply (3)

Description	Information
X016-1	Control panel connection
X018-1	Drive motor limit switches (1) - (3)
X019-1	Frequency converter controller
X020-1	Safety chain and relay control
X021-1	JK coin device and external coin device connection
X024-1	Drive motor (1) voltage supply
X025-1	Drive motor (2) voltage supply

Not shown:

2-pin BU13	Drive motor (1) voltage supply
4-pin BU23	Drive motor (1) controller
2-pin BU14	Drive motor (2) voltage supply
4-pin BU24	Drive motor (2) controller
2-pin BU15	Drive motor (3) voltage supply
4-pin BU25	Drive motor (3) controller

Battery CR 2032 3 V

13144 / 0

- Pull out the strip before commissioning.

8.3 Wellsystem MEDICAL / MEDICAL_PLUS: EMC information

8.3.1 Table 201: Guidance and manufacturer's declaration – electromagnetic emissions

The massage unit has been designed for use in the electromagnetic surroundings listed below. The operator of this equipment must ensure that it is only used in one of these surroundings.

Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The massage unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The massage unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-11	complies	Obligations according to Subclause 6.2.2 of EN 61000-3-11:2000 have to be considered. Especially the line impedance of 0.14 Ohms must be guaranteed for pulse massaging.

8.3.2 Table 202: Guidance and manufacturer's declaration – electromagnetic immunity

The massage unit has been designed for use in the electromagnetic surroundings listed below. The operator of this equipment must ensure that it is only used in one of these surroundings.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact	±6 kV contact	–
	±8 kV air	±8 kV air	
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	–
Surge 61000-4-5	±1 kV differential mode	±1 kV differential mode	–
	±2 kV common mode	±2 kV common mode	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	–
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle	<5 % U_T (>95 % dip in U_T) for 0.5 cycle	–
	40 % U_T (60 % dip in U_T) for 5 cycles	40 % U_T (60 % dip in U_T) for 5 cycles	
	70 % U_T (30 % dip in U_T) for 25 cycles	70 % U_T (30 % dip in U_T) for 25 cycles	
	<5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 5 s	

Note: U_T is the a.c. mains voltage prior to application of the test level.

8.3.3 Table 204: Guidance and manufacturer's declaration – electromagnetic immunity

The massage unit has been designed for use in the electromagnetic surroundings listed below. The operator of this equipment must ensure that it is only used in one of these surroundings.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Guided HF interference IEC 61000-4-6	3 Vrms 150 kHz up to 80 MHz	3 Vrms	<p>Portable and mobile wireless equipment is not to be used closer to the massage unit, including the lines, than the recommended safety distance, which has been calculated for the transmitting frequency using a suitable equation.</p> <p>Recommended safety distance:</p> $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz up to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz up to } 2.5 \text{ MHz}$
Radiated HF interference IEC 61000-4-3	3 V/m 80 MHz up to 2.5 GHz	3 V/m	<p>Where P is the transmitter's maximum output power (rated power) in Watts (W) in accordance with the details supplied by the manufacturer of the transmitter and the recommended safety distance is given in metres (m).</p> <p>The field strength of the stationary radio transmitter is less than the agreed level for all frequencies^b as determined by an on-site test^a.</p> <p>The following symbol must be used in the areas surrounding the equipment to indicate that interference is possible.</p> 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the massaging unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the unit.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

8.3.4 Table 206: Recommended separation distances between portable and mobile RF communications equipment and the massage unit

The massage unit has been designed for use in electromagnetic surroundings, in which the radiated HF interference is controlled. The operator of the massage unit can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF communications devices (transmitters) and the unit, in accordance with the recommended maximum output power of the communications devices as shown below.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

9 Index

A

Aligning the massage unit	
Empty well	36
Assembling the frame	30, 46

C

Coin devices	59
Commissioning	12
Connecting up the cooling system	45
Connecting up the pump	33
Connection ratings	58
Control box cable	34
Controller	61
Cooling system	45
Customer service	2

D

Dimensions	
MEDICAL	56
MEDICAL_PLUS	57
RELAX	54
RELAX_PLUS	55
Disassembling the frame	28
Disconnect the device from the power supply and ensure that it cannot be switched on again	13
Disposal	14

E

Electrical connections	11
EMV information (MEDICAL / MEDICAL_PLUS)	63
Equipment	15


F

Filling the well	43
Fitting the control box	34
Fitting the control panel	
MEDICAL / MEDICAL_PLUS	49
RELAX / RELAX_PLUS	48
Fitting the drive motors	44
Fitting the handle	47
Fitting the panels	51
Fitting the well	31

I

Information for service staff	8
Installation location	9

M	
Maintenance.....	12
N	
Nozzle carriages	
Check whether the nozzle carriages move smoothly	37
P	
Packaging	14
Packaging units.....	16
Power	58
R	
Recycling.....	14
Removing the control unit.....	21
Removing the drive motors	24
Removing the frame.....	17
Removing the panels.....	17
Removing the pump	25
Removing the well.....	25
S	
Safety information for dismantling and assembly.....	6
Scope of delivery.....	15
Stickers fitted on the device	60
Symbols	5
T	
Table of contents.....	3
Technical data.....	54
Test run	12
Time control	
accessories	59
Tools	15
Transport and transport damage.....	8
W	
Warning sticker	60
Water supply	12
Water/cooling conditions	59
Weights	58

1009892

Wellsystem GmbH

Koehlershohner Strasse
53578 Windhagen/ Germany
Tel. 02224/818-257
Fax 02224/818-254
info@wellsystem.de
www.wellsystem.com

Wellsystem GmbH

Koehlershohner Strasse
53578 Windhagen/ Germany
Phone +49 2224/818-250
Fax +49 2224/818-254
contact@wellsystem.com
www.wellsystem.com

wellsystem™