

Hochwertiger Naturkautschuk


Für die weiche Liegefläche unserer Wellsystem Wave_Touch und Medwave_Touch verwenden wir ausschließlich Naturkautschuk in Kombination mit einer hochfesten und besonders pflegeleichten Oberfläche. Wie bei allen Naturprodukten sind kleine Unebenheiten aufgrund des Produktionsprozesses unvermeidbar. Diese stellen jedoch keinen Reklamationsgrund dar. Keine Sorge: Die mechanischen Eigenschaften und die Haltbarkeit der weichen Liegefläche werden dadurch nicht beeinträchtigt.



Kleine Mulden



Kleine Blasen



Glatte Stellen

High-quality natural rubber


For the soft lying surface of our Wellssystem Wave_Touch and Medwave_Touch, we exclusively use natural rubber combined with a high-strength and particularly easy-care surface. As with all natural products, small irregularities due to the production process are unavoidable. However, these do not constitute a reason for complaint. No worries: The mechanical properties and durability of the soft lying surface are not affected by this.



Small indentations



Small bubbles



Smooth spots